

GEISTER

GEISTER®

The better way to operate™




Articulated
Instruments
for MIS




10-0255 Titanium Resano Dissecting Forceps 28cm

Optimised grasping surface with parallel grasping capability for a mini-thoracotomy approach this instrument helps to operate in common manner through smaller incisions. It provides strong but less traumatic grasping for valves and other vascular structures, as the pressure of jaws is equally distributed on the tissues. Double-action jaws.




26-4000 Titanium Papillary Muscle Scissors 28cm straight
26-4001 Titanium Papillary Muscle Scissors 28cm 85° angled

Specially designed for use in the papillary muscle plastics with extra wide jaw opening. Double-action jaws.



24-0540 Titanium De Bakey Needle Holder 28cm

Provides maximum needle grasping capability with Sapphire™ jaws for a mini-thoracotomy approach this instrument helps to operate in common manner through smaller incisions and reduces fatigue during long procedures. Recommended for needle sizes 4/0 to 6/0. Double-action jaws.



© Geister 2009 · Printed in Germany · All rights reserved.



The better way to operate™

Geister Medizintechnik GmbH
Föhrenstrasse 2
D-78532 Tuttlingen/Germany
Phone +49 (0) 7461 – 966 240
Telefax +49 (0) 7461 – 966 2422
E-mail info@geister.com
Internet www.geister.com

Distributor