


GEISTER

GEISTER®

The better way to operate™



Jacobson Soft Touch
Micro Needle Holders

Microsurgery redefined:



that special soft touch

- A new softer ratchet
- 40% increased needle grip
- perfect balance
- better control and feeling
- less fatigue
- faster action and response
- precise anastomosis

New GEISTER Jacobson Soft Touch Micro Needle Holders with newly designed softer ratchets and springs. Ratchets are easy to close and to release with minimum force and provide faster action and response for better control during anastomosis and less fatigue for hand muscles. Precise anastomosis results can be achieved.

Comes with Sapphire™ enhanced needle grip surface on its jaws for most secure needle grip.

Order information:

03-0500	1.2 mm tip, 18 cm, straight, with catch
03-0501	0.8 mm tip, 18 cm, straight, with catch
03-0502	0.4 mm tip, 18 cm, straight, with catch
03-0505	1.2 mm tip, 18 cm, curved, with catch
03-0506	0.8 mm tip, 18 cm, curved, with catch
03-0507	0.4 mm tip, 18 cm, curved, with catch
03-0510	1.2 mm tip, 21 cm, straight, with catch
03-0511	0.8 mm tip, 21 cm, straight, with catch
03-0512	0.4 mm tip, 21 cm, straight, with catch
03-0515	1.2 mm tip, 21 cm, curved, with catch
03-0516	0.8 mm tip, 21 cm, curved, with catch
03-0517	0.4 mm tip, 21 cm, curved, with catch

1.2 mm tip recommended for needles 5-0 and smaller,
0.8 mm tip recommended for needles 6-0 and smaller,
4.0 mm tip recommended for needles 8-0 and smaller.
Optionally available ceramic coated with „Black Diamond Film“.
Versions without catch on inquiry.



GEISTER®

The better way to operate™

Geister Medizintechnik GmbH
Foehrenstraße 2
D-78532 Tuttlingen/Germany
Phone +49-(0)74 61-96 62 40
Telefax +49-(0)74 61-96 62 422
E-mail: info@geister.com
Internet: www.geister.com

Distributor