

GEISTER



GEISTER®

The better way to operate™

Dingo-type
Bulldog Clamps



These reusable, featherweight soft bulldog clamps are especially suitable for temporary vessel occlusion in beating-heart surgery. Damage and drag to skeletonised arterial grafts can be reduced due to their softness and low weight.

The anti-magnetic clamps are made from titanium (grade 5) and are colored by anodization. Unlike steel, titanium does not stiffen with use, so these carefully calibrated titanium arterial clamps shall remain so indefinitely. They can be applied and removed with a utility forceps.



Order information:

Art.-No.	Overall length	Jaw width	for vessels diameter
20-0220	37.5mm	2.0mm	4mm
20-0221	30.0mm	2.0mm	4mm
20-0222	24.0mm	1.5mm	3mm
20-0223	14.0mm	1.0mm	2mm
20-0224	11.5mm	0.6mm	1mm



The better way to operate™

Geister Medizintechnik GmbH
Foehrenstraße 2
D-78532 Tuttlingen/Germany
Phone +49-(0)74 61-96 62 40
Telefax +49-(0)74 61-96 62 422
E-mail: info@geister.com
Internet: www.geister.com

