

GEISTER



GEISTER®

The better way to operate™

Xenon Powered
Headlight
Systems



Turn on the light !

Xenon Powered Headlight



Whisper quiet fan,
low heat development

Xenon powered with
daylight quality 5800 K

Movable stand with
dual hangers

Lightweight headlight system,
with dual adjustment headband,
neoprene fabric comfort pads,
bifurcated fiberoptic cable,
universal screw connector with
Storz-type adapter included,
sterilizable joystick, extended
linkage and variable spot size.

Order information:

- 02-1020 Bifurcated Headlight System with universal connector and Storz-type lightsource adapter, joystick, 2 gown clips, in carrying case
- 33-0311 Wolf-type lightsource adapter
- 33-0312 ACMI-type lightsource adapter
- 33-0313 Olympus-type lightsource adapter
- 02-1100 Movable Stand with hangers
- 33-0130 LISA smart Xenon Illuminator 100W 5800K with Storz-type socket (other sockets on inquiry)
- 33-0121 Replacement Bulb 100W Xenon (lifetime approx. 500 hours)
- 02-1032 Replacement Fiberoptic Cable 10 feet
- 02-1034 Replacement Joystick, sterilizable
- 02-1035 Replacement Sweat Pad, neoprene
- 02-1045 Replacement Gown Clips
- 02-0900 Mountable Multimedia-Kit CCD Camera System PAL
- 02-0901 Mountable Multimedia-Kit CCD Camera System NTSC



REF 02-1020



REF 33-0130

REF 02-1100



REF 02-0900

REF 02-1020 with 02-0900



Dual knob adjustment
headband

Moisture absorbing neoprene
fabric comfort pads

Sterilizable joystick and
extended linkage

Variable spot size
20-100 mm

Headband mountable
Camera System available
optionally

Universal screw connector,
Storz-type adapter included
(others optional)

Bifurcated 10 feet
transparent
cable



The better way to operate™



Light and comfort

Variable spot size
20-100 mm

Moisture absorbing
neoprene fabric

Softband for maximum
comfort

Order information:

- 02-1050 Bifurcated Comfort Softband Headlight System with universal connector and Storz-type lightsource adapter, joystick, 2 gown clips, in carrying bag
- 33-0311 Wolf-type lightsource adapter
- 33-0312 ACMI-type lightsource adapter
- 33-0313 Olympus-type lightsource adapter
- 02-1100 Movable Stand with hangers
- 33-0130 LISA smart Xenon Illuminator 100W 5800K with Storz-type socket (other sockets on inquiry)
- 33-0121 Replacement Bulb 100W Xenon (lifetime approx. 500 hours)
- 02-1032 Replacement Fiberoptic Cable 10 feet
- 02-1034 Replacement Joystick, sterilizable
- 02-1055 Replacement Softband, neoprene
- 02-1045 Replacement Gown Clips
- 02-0900 Mountable Multimedia-Kit CCD Camera System PAL
- 02-0901 Mountable Multimedia-Kit CCD Camera System NTSC



REF 02-1050



REF 33-0130

Lightweight softband headlight system, with neoprene fabric comfort headband, bifurcated fiberoptic cable, universal screw connector with Storz-type adapter included, sterilizable joystick, extended linkage and variable spot size.

GEISTER®
The better way to operate

Geister Medizintechnik GmbH
Foehrenstraße 2
D-78532 Tuttlingen/Germany
Phone +49-(0)74 61-96 62 40
Telefax +49-(0)74 61-96 62 422
E-mail: info@geister.com
Internet: www.geister.com

Distributor