

GEISTER®

The better way to operate™



Accu
Sternum
Saw

GEISTER

... the power of engineering

Accu Sternum Saw ... 9,6V NiMh-Power Technology



REF 14-6800

Accu Sternum Saw complete, with solid aluminium housing, brushless motor with 14000/min, 150W output, speed continuously adjustable, sensor control, includes 2 batteries, 2 saw blades, adjustable safety guard, battery charger, sterile hopper, carrying case.



REF 14-6804 Battery Ni-MH 9,6V 2700mAh
(approx. min. 250 times rechargeable)



REF 14-6802
Saw Blade for Sternum Saw

REF 14-6803
Safety Guard for Sternum Saw

REF 14-6805 Battery Charger 15V 18VA

GEISTER[®]
The better way to operate™

Geister Medizintechnik GmbH
Föhrenstraße 2
D-78532 Tuttlingen/Germany
Telefon +49 (0) 74 61 - 96 62 40
Telefax +49 (0) 74 61 - 96 62 422
E-mail info@geister.com
Internet www.geister.com

Distributor